## BIOTEX

(POLITETRAFLUOROETILENE)

Sutura sintetica non assorbibile

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### BIOTEX (POLYTETRAFLUROETHYLENE)

Synthetic non-absorbable suture

**DESCRIPTION**
BIOTEX sutures are made of monofilament non-absorbable polytetrafluoroethylene, a material with a porous microstructure, of part of the volume is composed of air.
BIOTEX sutures are available not coloured.
BIOTEX is available in a range of gauge sizes and lengths, non-needled or attached to stainless steel needles of different types and sizes. The needles' surface is provided with a silicone coat to make the tissues penetration easier and less traumatic for the patient.
BIOTEX sutures are available with pledgets in PTFE.
BIOTEX complies with the requirements of the European Pharmacopoeia's monograph "Non-absorbable sterile sutures".
BIOTEX is natural rubber latex free.

**Indications**
BIOTEX is indicated in all soft tissue approximation applications (suture with pre-assembled needle and) or/ ligature (suture with pre-assembled needle and free threads), in cardiovascular surgery, oral surgery, general surgery and in dura mater repair.

It is recommended for use where reduced suture line bleeding and to reduce hemostasis times, thanks to the porous microstructure of the thread which allows a quick attachment of the tissues and thanks to a 1:1 ratio between needle and thread.

Moreover, thanks to the smooth surface of the thread, BIOTEX allows less friction with the fabrics and greater maneuverability and smoothness, allowing precise positioning of the knots, an important feature in cardiovascular surgery and oral surgery.
BIOTEX is indicated, also, for the repair or replacement of chordae tendinae.

**Applications**
BIOTEX should be selected and implanted depending on patient conditions, surgical experience, surgical technique and wound size.

**Properties**
BIOTEX is not absorbable and it is chemically inert; the endo proliferation of the tissue in the suture thread may cause the suture thread to adhere to the tissue. This adhesion can make it difficult to remove the BIOTEX suture thread.

**Contraindications**
The use of BIOTEX is contraindicated in ophthalmic surgery, microsurgery and peripheral neural tissue.
The use of free threads is under no circumstances recommended in combination with other devices.

**Warnings/Precautions/Interactions**
The safety and effectiveness of this suture in peripheral neural, microsurgical and ophthalmic applications has not been established. The endo proliferation of the tissue in the suture thread may cause the suture thread to adhere to the tissue. This adhesion can make it difficult to remove the BIOTEX suture thread.

BIOTEX sutures must be exclusively used by specialized users, qualified for the surgical practice.
The correct application procedure reduces the risk of damages. The surgeon has the responsibility to select the correct format to use in different situations in order to minimize the risk of damages.
BIOTEX's available commercial variants can meet different needs. Before using BIOTEX, the user must verify the integrity of the package.
Open the primary package just a moment before the use, respecting asepsis' rules. While removing the package, then, the operator must make sure to preserve the sterility of the product contained in it. Take out the inner sterile package while using sterile gloves. The operator must have a good knowledge of surgical procedures in which BIOTEX is employed.

BIOTEX should be used, following the directions, in accordance with the usual suture and knotting's techniques and on the basis of the user's experience.

The risk of wound's dehiscence varies depending on the location of the wound and the suture employed. For this reason, the user must be experienced in surgical techniques that imply the use of Polytetrafluoroethylene sutures.
When using a suture, pay attention not to damage it, avoiding to crush or crimp it with surgical instruments such as forceps or needle-holders. Pull the wire between the needle's penetration channel and the needle itself. Don't pull the wire with too much strength. Avoid friction with gloves to prevent wire's surface abrasion. Pay attention not to cause damages and not to get hurt while handling the needle. The needle has to be grabbed between 1/3 and 1/2 of the distance from the attachment end to the point. If you grab it too close to the needle's eye, you can interfere with its penetration ability and it could break. If you take it too close to the

needle's eye, instead, it could bend or break. Trying to straighten up a bent needle, may induce its loss in bending and breaking resistance. Use needles must be treated following the safety procedures for the disposal of hospital wastes.

**Side effects**
Side effects associated with this product, are represented by an initial tissue inflammatory reaction of small entity and a temporary local irritation of the wound. As any foreign body, BIOTEX suture can aggravate a pre-existent wound.

**Sterility**
BIOTEX is sterilized with Ethylene Oxide. Do not re-sterilize. Do not use if the package is open or damaged.

**Storage**
It is recommended to store the device in a dry place, away from light sources. Do not use after the expiration date.

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