BIOCRYL
(POLYGLACTIN 910)
SUTURA STERILE SINTETICA ASSORBIBILE

DESCRIZIONE

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DESCRIPTION

Absorbable, sterile and synthetic BIOCRYL surgical sutures are made of a copolymer composed by glycolide for 90% and L- Lactide for 10%. They are braided and multifilament, coated by a mixture of calcium stearate and a copolymer of glycolide (30%) and lactide (70%). Copolymers used to manufacture this product are non-antigenic, non-pyrogenic and elicit only a slight tissue reaction during the absorption process. BIOCRYL sutures are natural rubber latex free.

BIOCRYL sutures are available both colored with Violet D&C N.2 dye (C.I. n.60725) and undyed. Moreover, they are available in a range of gauge sizes and lengths, non-needled or attached to stainless steel needles of different types and sizes. Needles' surface is provided with a silicone coat to make the tissues' penetration easier and less traumatic for the patient.

BIOCRYL complies with the European Pharmacopoeia monograph's requirements for absorbable sterile sutures.

INDICATIONS

BIOCRYL sutures are intended for use in general soft tissues approximation and/or ligation, including use in ophthalmic surgery, peripheral nerve anastomosis and microsurgery for vessels less than 2 mm diameter. The safety and effectiveness of BIOCRYL sutures in cardiovascular tissues have not been established.

APPLICATIONS

Sutures should be selected and implanted depending on patient conditions, surgical experience, surgical technique and wound size.

PROPERTIES

BIOCRYL sutures elicit a minimal initial inflammatory reaction in tissues and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of BIOCRYL sutures occurs by means of hydrolysis, where the copolymer degrades to lactic acids and glycolic acid which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength, followed by a loss of mass. Absorption of BIOCRYL sutures is approximately complete between 56 and 70 days. Studies on implantations on rats (performed with USP 2-0 and USP 4-0 sutures) showed that BIOCRYL sutures kept approximately the 75% of the original breaking strength after 2 weeks from the implantation. After three weeks, about 50% of the breaking strength is maintained. After 4 weeks the breaking strength was approximately the 30% of the original one for the USP 4-0 sutures and less than 25% for the USP 2-0 sutures.

Days implantation	Approximate percentage original strength remaining
14 days	75%
21 giorni	50%
28 days (USP 2-0)	23%
28 days (USP 4-0)	30%

CONTRAINDICATIONS

BIOCRYL sutures, being absorbable, should not be used where extended approximation of tissues under stress is required.

WARNINGS/ PRECAUTIONS/ INTERACTIONS

BIOCRYL sutures must be exclusively used by specialized users, qualified for the surgical practice.

Before using BIOCRYL, the operator must verify the integrity of the package.

Open the primary package just a moment before the use, respecting asepsis' rules.

While removing the wrapping, then, the operator must preserve the sterility of the product contained in it. Take out the inner package while wearing sterile gloves.

As any suture's material, pay attention not to damage the wire during handling, in particular it should not be crimped or crushed with forceps or needle-holders. While tightening the suture, pull the wire between the needle's penetration channel and the needle itself. Don't pull the wire with too much strength or over sharp objects. Avoid friction with gloves to prevent the wire's surface abrasion. The needle has to be grabbed at a point situated between 1/3 and 1/2 of the distance from the attachment end to the point. If you take it too near to the tip, you can interfere with its penetration ability e it could break. If you grab it too near to the needle's eye, instead, it could bend or break. Pay attention while using it in order to prevent the risk of needle stick injury. Used needles must be disposed according to the safety procedures for the disposal of hospital wastes.

The risk of dehiscence varies depending on the sites of application and the suture employed. The user has to be experienced in surgical procedures in which BIOCRYL is used. BIOCRYL, being an absorbable suture, can transiently act as a foreign body. As any foreign body, prolonged contact of the suture with salt solutions, such as those of urinary and biliary tracts, can induce the formation of calculi. Appropriate surgical procedures must be used for treating contaminated or infected wounds.

As this is an absorbable suture, the surgeon should eventually use supplementary non-absorbable sutures while closing sites that could expand, stretch or spread or request an additional support. Skin sutures that have to stay in place more than 7 days, may induce local irritations. For this reason, if necessary, the suture has to be removed. The use of absorbable sutures on tissues with a low vascularization should be carefully evaluated, because extrusion or delayed absorption may occur. Subcutaneous sutures have to be implanted as deep as possible in order to minimize the formation of erythema and hardening which normally occur during the absorption process.

BIOCRYL is not recommended in elderly, malnourished and debilitated patients, or in patients suffering from conditions which may delay wound healing.

SIDE EFFECTS

Side effects associated to the use of this product are represented by an initial tissue inflammatory reaction of small entity and a temporary local irritation of the wound. As any foreign body, BIOCRYL suture may aggravate a pre-existent infection.



STERILITY

BIOCRYL is sterilized with Ethylene Oxide. Do not re-sterilize. Do not use if the package is open or damaged.

STORAGE

It is recommended to store the device in a dry place, away from light sources. Do not use after the expiration date.

LEGEND OF THE SYMBOLS ON THE PACKAGE

	The product has to be stored in a dry place
	Keep the product away from the light
	Do not use if the package is damaged
	Carefully read the instructions for use
	Do not re-use
	Do not re -sterilize
	Reference code
	Lot number
	Expiration date (use within the date indicated, month and year)
	Manufacturer (Bio Sud Medical Systems S.r.l.)
	Sterile product. Sterilization method: Ethylene Oxide
	Colored, multifilament, braided, coated, absorbable suture
	Undyed, multifilament, braided, coated, absorbable suture
	Needle's typology: curvature, size and shape of the tip (an example is shown)
	Suture size following the United States Pharmacopoeia (USP)
	and the European Pharmacopoeia (EP – Gauge number)
	Suture's length

Data-matrix, reporting the reference to the product code, the lot number and the expiry date of the medical device

CE marking and identification number of the notified Organism. The product complies with the essential requirements of the Medical Devices Directive 93/42/CEE

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- Producător (Bio Sud Medical Systems S.r.l.)